

K103686  
7/15

21<sup>st</sup> January 2011  
Page: 01 of 04

MAR 15 2011

**510(K) SUMMARY**  
(Per section 807.92 ©)

|                                      |                |  |                                |
|--------------------------------------|----------------|--|--------------------------------|
| <b>CONTACT DATA</b>                  |                |  |                                |
| <b>Submitter's Name</b>              |                | Larsen & Toubro Limited  |                                |
| <b>Address</b>                       |                | L&T Medical Equipment & Systems,<br>Mysore Campus, Gate No.5, Plot No. 358-360,<br>KIADB Industrial Area, Hebbal,<br>Mysore – 570018, Karnataka, INDIA |                                |
| <b>Telephone</b>                     | 91-821-2407200 | <b>Fax</b>   | 91-821-2407001                 |
| <b>Contact Person</b>                | A.B.Deshpande  | <b>Title</b>   | Head – QA & Regulatory Affairs |
| <b>E-Mail address</b>                |                | DeshpandeAB@myw.ltindia.com  |                                |
| <b>Date the summary was prepared</b> |                | 21 <sup>st</sup> January 2011  |                                |

K103686  
#2/5

21<sup>st</sup> January 2011

Page: 02 of 04

| <b>DEVICE</b>              |                           |
|----------------------------|---------------------------|
| <b>Trade name</b>          | STAR 50N                  |
| <b>Common name</b>         | Patient Monitoring System |
| <b>Classification name</b> | Vital Signs Monitor       |

| <b>PREDICATE DEVICE IDENTIFICATION</b>         |   |                                |     |
|--|---|--------------------------------|-----|
| <b>CFR21 Section</b>                           | 870.2300  | <b>Product code (optional)</b> | MWI |
| <b>Classification panel</b>                    | Cardiovascular  |                                |     |
| <b>Device Class</b>                            | Class II  |                                |     |
| <b>Legally marketed Comparison Device / K#</b> | <ul style="list-style-type: none"><li>• STAR 55 Patient Monitoring System (L&amp;T Medical Equipments &amp; systems) / K080173</li><li>• Stellar 300 Patient Monitoring System (L&amp;T Medical Equipments &amp; systems) / K093017</li></ul> |                                |     |

21<sup>st</sup> January 2011

### **DEVICE DESCRIPTION**

STAR 50N is a multi-parameter patient monitoring system for continuous monitoring of the physiological parameters ECG (3/5 lead), Respiration, NIBP, IBP, Temperature, SpO2 and external CO2 (optional).

STAR 50N is a 6-channel monitor with 10.4" TFT display capable of displaying ECG, Respiration, Spo2, CO2, digital values of HR/PR, SpO2, RR, Non-Invasive Blood pressure (Systolic, Diastolic and Mean), Invasive Blood Pressure (Systolic, Diastolic and Mean), Temperature, EtCO2 and FiCO2 readings. It has 168 Hours tabular and graphical trends. It has special feature of NIBP having a trend of storing last 240 readings. It has Alarm Recall facility with last 24 patient alarms details. It has a two-channel thermal array recorder for printing of Tabular trends & waveforms. The device permits patient monitoring with adjustable alarm limits as well as visible and audible alarm signals. It has got optional communication features - USB, Wi/Fi (optional) and Ethernet. This monitor can also be connected to L&T Central Nursing Station (Skyline 55) and an external LCD-TFT display

### **INTENDED USE OF THE DEVICE**

The STAR 50N multi-parameter Patient Monitoring system is intended to monitor a single adult, pediatric and neonatal patient's vital signs at the bedside or during intra-hospital transport along with the appropriate accessories mentioned / supplied with the unit. Vital signs parameters include ECG (3 lead /5 lead), SpO2, Respiration, Temperature, external optional Capnography (CO<sub>2</sub>). It can also display the digital values of HR/PR, SpO<sub>2</sub>, RR, Non-Invasive Blood Pressure (Systolic, Diastolic and Mean), Invasive Blood Pressure (Systolic, Diastolic and Mean), Temperature, EtCO<sub>2</sub>, and FiCO<sub>2</sub>. This monitor can also be connected to L&T Central Nursing Station (Skyline 55) and an external LCD-TFT display.

The user, responsible to interpret the monitored data made available, will be a professional health care provider. The device permits patient monitoring with adjustable alarm limits as well as visible and audible alarm signals. The monitor is not intended for home use.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

**TECHNOLOGICAL CHARACTERISTICS OF THE DEVICE COMPARED TO THE PREDICATE  
DEVICE**

**Predicate Devices:**

| Company Name        | Trade Name  | 510(k) Number |
|---------------------|-------------|---------------|
| Larsen & Toubro LTD | Star 55     | K080173       |
| Larsen & Toubro LTD | Stellar 300 | K093017       |

**Predicate device:**

Star 55 is a Multi parameter patient monitor used as a predicate device for all the parameters of Star 50N Multi parameter patient monitor, except for NIBP parameter.

Stellar 300 is a three parameter (NIBP, SpO2 & single channel Temperature) patient monitor used as a predicate device for 'NIBP parameter' of Star 50N Multi parameter patient monitor.

Comparison of all the parameters of STAR 50N to that of the predicate devices is given in the "Substantial Equivalence Equipment comparison" document.

**Summary of Critical Parameters**

| Parameter                     | Star 50N | Star 55 |
|-------------------------------|----------|---------|
| 510(k) Number                 |          | K080173 |
| Alarms                        | Same     | Same    |
| ECG                           | Same     | Same    |
| Respiration                   | Same     | Same    |
| Temperature                   | Same     | Same    |
| Pulse Oximetry                | Same     | Same    |
| IBP (Invasive Blood Pressure) | Same     | Same    |
| Capnography                   | Same     | Same    |

| Parameter     | Star 50N | Stellar 300 |
|---------------|----------|-------------|
| 510(k) Number |          | K093017     |
| NIBP          | Same     | Same        |

**Compliance to standards:**

The following international standards are referred.

IEC 60601-1 Medical Electrical safety

IEC 60601-1-2 EMC compliance

**Conclusion:**

Based on the Technological characteristics of STAR 50N and its comparison with that of predicate devices Star 55 and Stellar 300, Larsen & Toubro Limited believes that their device is substantially equivalent to this predicate Monitors and doesn't pose any additional risk on safety & effectiveness of the device.

(N Ravindran)  
Head – Design & Development



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room -WO66-G609  
Silver Spring, MD 20993-0002

Larsen & Turbo Limited  
c/o Mr. E. J. Smith  
Regulatory Consultant  
Smith Associates  
1468 Harwell Avenue  
Crofton, MD 21114

MAR 15 2011

Re: K103686  
Trade/Device Name: Star 50N  
Regulatory Number: 21 CFR 870.2300  
Regulation Name: Physiological Patient Monitor (without arrhythmia detection or alarms)  
Regulatory Class: II (two)  
Product Code: 74 MWI  
Dated: February 25, 2011  
Received: February 28, 2011

Dear Mr. Smith:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Page 2 – Mr. E. J. Smith

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



Bram D. Zuckerman, M.D.  
Director  
Division of Cardiovascular Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

# Indications for Use Form

## Indications for Use

510(k) Number (if known): K103686

Device Name: Star 50N

### Indications for Use:

The Star 50N multi-parameter Patient Monitoring system is intended to monitor a single adult, pediatric and neonatal patient's vital signs at the bedside or during intra-hospital transport along with the appropriate accessories mentioned / supplied with the unit. Vital signs parameters include ECG (3 lead / 5 lead), SpO2, Respiration, Temperature, external optional Capnography (CO2). It can also display the digital values of HR/PR, SPO2, RR, Non-Invasive Blood Pressure (Systolic, Diastolic and Mean), Invasive Blood Pressure (Systolic, Diastolic and Mean), Temperature, EyCO2, and FiCO2.

The user, responsible to interpret the monitored data made available, will be a professional health care provider. The device permits patient monitoring with adjustable alarm limits as well as visible and audible alarm signals. The monitor is not intended for home use.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Prescription Use   X   AND/OR Over-The-Counter Use \_\_\_\_\_  
(Part 21 CFR 801 Subpart D) (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

*[Handwritten Signature]*  
*[Handwritten Signature]*  
*3/15/2011*

Page   1   of   1  

(Division Sign-Off)  
Division of Cardiovascular Devices

510(k) Number   K103686